Technology and Health Care -1 (2023) 1–12

DOI 10.3233/THC-231006

IOS Press

# Graphic type differentiation of cell count data for diagnosis of early and late periprosthetic joint infection: A new method

- Bernd Fink<sup>a,b,\*</sup>, Marius Hoyka<sup>a</sup>, Benedikt Blersch<sup>a</sup>, Hannsjörg Baum<sup>c</sup> and
- 5 Florian Hubert Sax<sup>a</sup>
- <sup>a</sup>Department for Joint Replacement, Rheumatoid and General Orthopaedics, Orthopaedic Clinic
- Markgröningen, Markgröningen, Germany
- 8 Orthopaedic Department, University Hospital Hamburg-Eppendorf, Hamburg, Germany
- <sup>9</sup> Institute for Laboratory Medicine and Transfusion Medicine, RKH Regionale Kliniken Holding und
- 10 Services GmbH, Ludwigsburg, Germany
- 11 Received 31 July 2023
- 12 Accepted 21 September 2023
  - Abstract.
- BACKGROUND: Graphic type differentiation of cell count data of synovial aspirates is a new method for the diagnosis of early and late periprosthetic joint infection.
- OBJECTIVE: The aim of the study was to analyse if the same 6 LMNE-types can be differentiated in the new Yumizen H500 cell counter as it was the case for the old cell counter ABX Pentra XL 80 of previous publications, to verify if the erythrocyte and thrombocyte curves of the new device give additional information and to calculate the difference of cell count in LMNE-type I and III (with abrasion) in the cell counter and in the manual counting chamber (Neubauer improved).
- METHODS: 450 aspirates of 152 total hip arthroplasties and 298 knee arthroplasties obtained for the diagnosis of periprosthetic joint infection were analysed with the Yumizen H500.
- RESULTS: All LMNE-matrices of the 450 aspirates could assigned to one of the six LMNE-types. There were 76 LMNE-type I, 72 LMNE-type II, 14 LMNE-type III, 241 LMNE-type IV, 36 LMNE-type V and 12 LMNE-type VI. The erythrocyte and thrombocyte distribution curves were very helpful for differentiation of hematoma and infection. The cell count in the manual counting procedure was lower than in the cell counter: for the LMNE-type I (abrasion type) the median of the difference was 925/µL (median) and for the LMNE-type III (combined type of infection and abrasion) 3570/µL (median).
- 27 **CONCLUSION:** The described graphic type differentiation is a new and helpful method for differentiation of hematoma and early PJI as well as abrasion and late PJI.
- Keywords: Leucocyte count, synovia, periprosthetic joint infection, cell counter

# 1. Introduction

Periprosthetic joint infection (PJI) is a devastating complication of arthroplasty procedures and has

<sup>\*</sup>Corresponding author: Bernd Fink, Department for Joint Replacement, Rheumatoid and General Orthopaedics, Orthopaedic Clinic Markgröningen, Kurt-Lindemann-Weg 10, 71706 Markgröningen, Germany. Tel.: +49 7145 9153201; E-mail: bernd.fink@rkh-gesundheit.de.

B. Fink et al. / Graphic type differentiation of cell count data for diagnosis of early and late PJI.

many consequences. The level of incidence for primary total hip and knee arthroplasties ranges between 1% and 2% on average [1–4]. With regard to revision arthroplasty, it is even higher and ranges between 2% and 7.5% [5–7]. The accuracy of the preoperative diagnosis of possible PJI is critical because, if it remains unrecognized before or during revision, it leads to a high risk of failure due to recurrence of the infection [8,9].

There is no single diagnostic test with an accuracy of 100% to verify or rule out an early/acute or late/chronic periprosthetic joint infection; thus, combinations of different diagnostic tests are helpful [10–13]. One important diagnostic method for PJI is the determination of the leukocyte count (WBC) in the joint synovia. Some authors consider it one of the most important diagnostic parameters [14–16] and it is included as one of the minor criteria in the definition of the PJI, the MSIS criteria and the more recent ICM criteria [10,17,18].

To set the cut-off value of WBC for early/acute and late/chronic periprosthetic joint infection is difficult for several reasons and therefore varies between the different published reports. In a consensus meeting, the threshold of WBC in the aspirate for acute periprosthetic infection (within 4 weeks after implantation or onset of symptoms in acute hematogenous infections) was set at 10,000 cells/μL, because different publications quoted threshold values between 8,910 cells/μL und 12,800 cells/μL [19–24]. However, this value corresponds to leukocyte numbers that may also be present in normal blood, so that postoperative hematomas in the joint may also exhibit leukocyte counts at this level [19,21]. Another important reported parameter is the percentage of polymorphonuclear leukocytes in the aspirate [10,17,18,23], whereby the cut-off value was determined to be 89% in several studies [19–21] and 79.5% in another study [24]. The International Consensus Group meeting fixed the cut-off value at 90% PMN in the aspirate [23]. However, the percentage of polymorphonuclear leukocytes in the aspirate can also be influenced by the presence of a hematoma. Because postoperative edema and swelling can make it difficult to distinguish clinically between early periprosthetic infection and hematoma, additional tools or parameters for the differentiation of these two situations would be helpful.

In chronic/late periprosthetic joint infections, the variance of thresholds for leukocyte number in the synovial aspirate is even higher [25]. Thresholds between 1,100 cells/ $\mu$ L [26] and 5,000 cells/ $\mu$ L [27] are described in the literature [16,25]. This leads to the setting of the threshold at > 3,000 cells/ $\mu$ L (where PJI is confirmed) [23,28], and < 1,500 cells/ $\mu$ L (where PJI is unlikely) [28]. Factors that lead to this broad variance of thresholds include the time after the operation, the duration of the symptoms, the detected microorganisms, administered antibiotics, and the authors only setting the cut-off value after receiving the results, depending on the prevalence of high sensitivity or specificity [29].

In a previous paper, we could show that this variance was also caused by abraded wear particles from the articulating surface of the artificial joint (polyethylene particles or metal particles) that were detected and counted by the cell counter in addition to the leucocytes in the aspirate [25]. We could differentiate between leucocytes and wear particles in the aspirate because of the different behaviour with regard to light absorption and we could also identify different graphical types (so called LMNE-types) of the aspirate: type I (abrasion type), type II (infection type), type III (mixed type of infection and abrasion), type IV (indifference type) [25]. We could show that these graphical LMNE types had a significantly high correlation to the histological types I–IV (with the same names) according to Krenn et al. [30,31] and Morawietz [32].

In a second paper, we expanded the number of graphical LMNE-types to include the differentiation of periprosthetic infection and hematoma: type V (hematoma) and type VI (mixed type of infection and hematoma) [33]. Differentiation of LMNE types into infection (types II and VI) and non-infection (types IV and V) resulted in a sensitivity of 100%, a specificity of 97.3%, and a positive likelihood ratio of 37.0 [33].




Fig. 1. a: LMNE matrix with the different areas corresponding to the leukocyte populations and the NOISE area of the ABX Pentra XL 80. b: LMNE matrix with the different areas corresponding to the leukocyte populations and the NOISE area of the Yumizen H500.

The cell counter ABX Pentra XL 80 (Horiba Medical, Montpellier, France) was used for the investigations described in those two papers. In the meantime, the next generation of these cell counters is on the market – the Yumizen H500 (Horiba Medical, Montpellier, France) – and has replaced the old ones. As a result, the LMNE-plots look somewhat different to those of the ABX Pentra XL 80 (Fig. 1a and b).

Therefore, the questions of the current study were:

- 1. Can the same 6 LMNE-types be differentiated in the new Yumizen H500 as was the case for the old ABX Pentra XL 80?
- 2. Does the new cell counter give additional information for the differentiation of these 6 types?
- 3. What difference of cell count are measured in LMNE-type I and III (with abrasion) in the cell counter and in the manual counting chamber (Neubauer improved).

# 2. Material and methods

78

79

80

81

82

83

84

85

86

87

88

89

450 aspirates of 152 total hip arthroplasties and 298 knee arthroplasties obtained for the diagnosis of periprosthetic joint infection were analysed with the Yumizen H500. There were 239 female and 211






Fig. 2. LMNE matrix of a type I with metal debris particles of a 84-year old female patient with knee arthroplasty. The cluster in the NOISE-area. The measured "cell count" was 1970 cells/ $\mu$ L and the percentage of polymorphnuclear neutrophil leukocytes (PMN) was 23.4%.

male patients with a mean age of  $70.5 \pm 16.8$  (range 50–90) years. Systemic inflammatory diseases like rheumatoid arthritis were excluded, because these diseases can have leukozytes in the joint beside a PJI [34]. The Yumizen H500 (Horiba Medical, Montpellier, France) is a laboratory diagnostic device for the analysis of the cell count and the WBC-differentiation of blood and also of other body fluids. Here we selected the so-called 5-DIFF mode from the various processing modes available. A total of 26 laboratory parameters are recorded, including the 5 cell types – eosinophils, neutrophils, monocytes, lymphocytes and basophils – as well as atypical lymphocytes and large, immature cells. These cell types are graphically mapped in a so-called LMNE matrix depending on their cell volume (x-axis) and their light scattering or refraction and absorption (y-axis) (Fig. 1b). The analysis is based on an impedance measurement, flow cytometry and cytochemistry. The volume difference of the cells due to the impedance measurement is shown graphically on the X-axis and the differentiation of the light absorption in flow cytometry on the Y-axis of the LMNE matrix. This enables the graphical assignment and thus differentiation of the 4 leukocyte populations: lymphocytes, monocytes, neutrophils and eosinophils (Fig. 1b). Impurities appear in the so-called NOISE area of the LMNE matrix; in our synovial analyses of joints with prostheses these are abraded wear particles (Fig. 1b). The LMNE-matrices of each patient were allocated to one of the different




Fig. 3. LMNE matrix of a type II (infection type) with a cluster of data points in the area of the neutrophil leukocytes and no peak in the erythrocyte field in a 57-year old patient with a late periprosthetic joint infection of a total knee arthroplasty. The measured cell count was 9840 cells/µL and 89.4% PMN.

LMNE-types (I to VI) (if possible). Additionally, the cell count and percentage of polymorphonuclear leucocytes were also measured by the Yumizen H500 cell counter. When a LMNE-type with abrasion (types I and III) was identified, the cell count was also determined manually using a Neubauer improved counter chamber (Carl Roth GmbH, Karlsruhe, Germany). Additionally the cell counter Yuminzen H500 gave curves for the erythrocyte and thrombocyte distribution of the aspirate with the size of the cell on the x-axis and the amount of cells on the y-axis (Figs 2–7).

The evaluations and assignment of the matrices to the different types were carried out twice by 2 examiners (BF and MH) independently of one another. IBM SPSS Version 24 (IBM Corp., Aemonk, NY, USA) and Microsoft Excel (Microsoft, Redmond, WA, USA) were used for statistical analysis. Categorical variables are depicted as frequencies, while continuous variables are shown as medians, standard deviation and ranges.

## 3. Results

107

108

109

110

111

112

113

114

115

116

117

118

All LMNE-matrices of the 450 aspirates could assigned to one of the six LMNE-types. There were

121

122

123

B. Fink et al. / Graphic type differentiation of cell count data for diagnosis of early and late PJI

### Table

Assignment of the aspirates to the different LMNE-types and additional identification of abrasion debris (substantial in type I and III and slight in the other types) as well as hematoma (substantial in type V and VI and slight in the other types)

LMNE-type	N	With abrasion	With hematoma
I (abrasion type)	76	76	3
II (infection type)	72	9	20
III (mixed type = infection $+$ abrasion)	14	14	2
IV (indifference type)	241	88	37
V (hematoma type)	36	7	36
VI  (mixed type = infection + hematoma)	12	3	12




Fig. 4. LMNE matrix of a type III (combined type of infection and abrasion) with a cluster in the area of the neutrophil leukocytes and a second cluster in the NOISE area as well as no peak in the erythrocyte field in a 83-year old female patient with a periprosthetic joint infection of a hip arthroplasty. The measured cell count was 2940 cells/µL with 79.6% PMN.

76 LMNE-type I, 72 LMNE-type II, 14 LMNE-type III, 241 LMNE-type IV, 36 LMNE-type V and 12 LMNE-type VI (Table 1). The assignment by the two examiners (BF, MH) showed a high reliability, with an intrarater intraclass correlation coefficient of 0.99 and of 0.98 between raters, respectively.

For the detection of hematoma, the erythrocyte and thrombocyte distribution curves were very helpful

125

126

127

128

129

Table 2
Cell counts and percentages of polymorphonuclear leucocytes
of the different types

I MNE trino	Cell count ( $/\mu$ L)	% PMN
LMNE-type	(median; min-max)	(median; min-max)
Type I	1,195; 274–158,070	24.6; 4.0–59.2
Type II	6,595; 1,130–179,510	81.2; 69.1–94.1
Type III	9,105; 1,620–93,290	69.5; 39.3–89.8
Type IV	560; 90–3,120	31.8; 3.7–77.9
Type V	2,740; 220–12,980	58.2; 20.6-89.0
Type VI	56,880; 6,380–287,630	80.8; 52.9–91.5




Fig. 5. LMNE matrix of a type IV (indeterminate type) with no obvious data point clusters or increase in cell types or particles in a 62-year old female patient. The measured cell count was 400 cells/ $\mu$ L with 22.5% PMN.

in indicating the existence of a significant number of erythrocytes and thrombocytes. This allowed the differentiation between the hematoma (type V) with a high number of leucocytes, and infection (type II) to be carried out more easily: the hematoma exhibited a peak in the erythrocyte and thrombocyte curve whereas the infection type did not (Figs 3 and 6). The higher the peak of the erythrocyte and thrombocyte curve, the more erythrocytes and thrombocytes (and more hematoma) were in the aspirate.

The cell counts and percentages of polymorphonuclear leucocytes of the different types are shown

Fig. 6. LMNE matrix of a type V (hematoma) 75-year old female patient with a hematoma 6 weeks after knee arthroplasty. There is a cluster in the field of the lymphocytes and another clear cluster in the field of the neutrophils and a peak in the erythrocyte field. The measured cell count" was  $10,120 \text{ cells}/\mu\text{L}$  with 59.2% PMN.

in Table 2. The cell count in the manual counting procedure was lower than in the cell counter: for the LMNE-type 1 (abrasion type) the median of the difference was  $925/\mu$ L with a range of 90 to 45,370 cells/ $\mu$ L and for the LMNE-type III (combined type of infection and abrasion)  $3570/\mu$ L (median) with a range of 1,200 to 10,210 cells/ $\mu$ L, respectively.

# 4. Discussion

130

131

132

133

134

135

136

137

138

139

140

141

142

The analysis of the aspirates by differentiating the LMNE-matrices into the LMNE-types I to VI with the ABX Pentra XL 80 (Horiba Medical, Montpellier, France), as described in the two previous publications, could clearly be reproduced with the new generation of this laboratory diagnostic device, the Yumizen H500 (Horiba Medical, Montpellier, France) [25,33]. Therefore, it can be expected that analysing joint aspirates with this new device will give additional information for the diagnostic of periprosthetic joint infections, as was the case for the older device.

On the basis of the results of an earlier paper, in late/chronic PJI, the additional use of the LMNE-type classification can help to distinguish between a true infection type and a debris type, which would




Fig. 7. LMNE matrix of a type VI (combined type II of hematoma and infection) 5 weeks after implantation of a hip arthroplasty with a cluster in the area of the neutrophil leukocytes and increases in the other areas of the white blood cells as well as a peak in the erythrocyte field in a 72 old female patient with an early periprosthetic joint infection. The measured cell count was 15200 cells/ $\mu$ L with 84.4% PMN.

otherwise lead to concluding an incorrectly high cell number [25]. The much lower number of cells counted in the manual counting chamber support the hypothesis that abraded wear particles can be mis-counted as leucocytes in the cell counter, especially in LMNE-type I and III aspirates. Above all, this method can be used to differentiate leukocytes from metallic abrasion particles, which appears to be of particular importance, since joint aspirates with metal abrasion can look like pus, can have very high cell counts, elevated CRP- and alpha-defensin values and can thus wrongly be interpreted as a periprosthetic infection [35–40]. Furthermore, the combination of cell count-analysis and LMNE types II or III-determination enabled us in the first paper to reduce the cut-off value for late PJI of the cell count in the aspirate to 1,400 cells/ $\mu$ L without losing high sensitivity (sensitivity of 90.2% and a specificity of 91.9%, receiver operation characteristics curve) [25]. The combination of cell counting and the graphic representation in the LMNE matrix means that fewer periprosthetic infections are overlooked and the diagnostic value of the cell count analysis in the joint aspirate is increased.

Based on the results of another previous paper, the LMNE matrix evaluation during the first weeks after implantation of an endoprosthesis can help to differentiate infections from hematomas, especially in the cases with cell counts in the borderline range of the threshold value chosen by the ICM. This

159

160

161

162

163

164

165

166

167

168

169

170

171

172

173

176

177

178

179

180

181

182

183

184

185

186

187

188

189

190

191

192

193

194

195

196

197

14/11/2023; 10:54

is because cluster plots in the LMNE-fields representing lymphocytes, basophils and/or eosinophils show that a significant part of the cluster in the neutrophil leukocyte field must be due to a hematoma. Therefore, cell count values in the borderline range can be corrected downward and thus not be interpreted as an infection. Hereby, the additional erythrocyte and thrombocyte distribution plot of the Yumizen H500 helps to identify a hematoma and distinguish it from an infection. Therefore, this method can be used to differentiate infections from hematoma. This is of particular importance because postoperative edema and swelling can make it difficult to distinguish these two situations clinically and the cell count in a hemarthrosis can be at the level of the threshold for infection described by the ICM and other authors [19–24].

This study has some limitations. First of all, only the reproducibility of the differentiation of the 6 LMNE-types with the new laboratory device was analyzed. The accuracy, sensitivity and specificity of the LMNE-type-differentiation with the Yumizen H500 have to be analyzed in further studies using criteria like that of the MSIS for comparison. Moreover, the type classification is somewhat dependent on the personal interpretation and experience of the examiner. Even if the reliability of this type classification was very high in our study, this does not rule out a certain subjectivity in the interpretation of the LMNE matrices.

# 5. Conclusion

The graphic representation of the cell count analysis of synovial aspirates of joints with endoprostheses is a new and helpful method for determining the cell count of the aspirate and to differentiate between true early periprosthetic infections with increased leukocyte count and hematomas in the early postoperative period. In addition, the ability to differentiate real late PJIs with increased leukocyte counts and apparently increased numbers due to particle abrasion increases the diagnostic value of the cell count analysis in such cases. In our opinion, it should therefore be included in the diagnostic armamentarium used to identify periprosthetic joint infections in future.

# Conflict of interest

None to report.

# References

- Li C, Renz N, Trampuz A. Management of periprosthetic joint infection. Hip Pelvis. 2018; 30: 138-146.
- Bohl DD, Sershon RA, Fillingham YA, Della Valle CJ. Incidence, risk factors, and sources of sepsis following total joint arthroplasty. J Arthroplasty. 2016; 31: 2875-2879.
- Kurtz SM, Lau E, Watson J, Schmier JK, Parvizi J. Economic burden of periprosthetic joint infection in the United States. J Arthroplasty. 2012; 27: 61-65.
- Pulido L, Ghanem E, Joshi A, Purtill JJ, Parvizi J. Periprosthetic join infection: The incidence, timing, and predisposing factors. Clin Orthop Relat Res. 2008; 466: 1710-1715.
- Fink B, Urbansky K, Schuster P. Mid-term results with the curved modular tapered, fluted titanium Revitan stem in revision hip replacement. Bone Joint J. 2014; 96-B: 889-95.
- Quinlan ND, Werner BC, Brown TE, Browne JA. Risk of prosthetic joint infection increases following early aseptic revision surgery of total hip and knee arthroplasty. J Arthroplasty. 2020; 35: 3661-3667.
- Zingg M, Kheir MM, Ziemba-Davis M, Menethini RM. Reduced infection rate after aseptic revision total knee arthroplasty with extended oral antibiotic protocol, J Arthroplasty, 2022; 37: 905-909.

220

221

222

223 224

225

226

227

228

229

234 235

236

240

241

247

248

249

250

251

252

253

- Della Valle CJ, Zuckermann JD, Di Cesare PE. Periprosthetic sepsis. Clin Orthop Rel Res. 2004; 420: 26-31. 198
  - Hanssen AD. Managing the infected knee: As good as it gets. J Arthroplasty. 2002; 17(4 Suppl 1): 98-101.
- 200 [10] Parvizi J, Tan TL, Goswami K, Higuera C, Della Valle C, Chen AF, Shohat N. The 2018 definition of periprosthetic hip and knee infection: An evidence-based and validated criteria. J Arthroplasty. 2018; 33: 1309-1314. 201
- Fink B, Schuster P, Braun R, Tagtalianidou E, Schlumberger M. The diagnostic value of routine preliminary biopsy in 202 diagnosing late prosthetic joint infection after hip and knee arthroplasty. Bone Joint J. 2020; 102-B: 329-335. 203
- Fink B, Makowiak C, Fuerst M, Berger I, Schäfer P, Frommelt L. The value of synovial biopsy, joint aspiration and 204 C-reactive protein in the diagnosis of late peri-prosthetic infection of total knee replacements. J Bone Joint Surg Br. 2008; 205 206
- Fink B, Gebhard A, Fuerst M, Berger I, Schäfer P. High diagnostic value of synovial biopsy in periprosthetic joint 207 208 infection of the hip. Clin Orthop Relat Res. 2013; 471: 956-964.
- Della Valle CJ, Sporer SM, Jacobs JJ, Berger RA, Rosenberg AG, Paprosky WG. Perioperative testing for sepsis before 209 revision total knee arthroplasty. J Arthroplasty. 2007; 22(6 Suppl 2): 90-93. 210
- Trampuz A, Hanssen AD, Osmon DR, Mandrekar J, Steckelberg JM, Patel R. Synovial fluid leukocyte count and 211 differential for the diagnosis of periprosthetic knee infection. Am J Med. 2004; 117: 556-562. 212
- [16] Jandl N, Kleiss S, Mussawy, Beil FT, Hubert J, Rolvien T. Absolute synovial polymorphonuclear neutrophil cell count as 213 a biomarker of periprosthetic joint infection. Bone Joint J. 2023; 105-B: 373: 381. 214
- Parvizi J, Zmistowski B, Berbari EF, Bauer TW, Springer BD, Della Valle CJ, Garvin KL, Mont MA, Wongworawat MD, 215 Zalavras CG. New definition for periprosthetic joint infection: From the Workgroup of the Musculoskeletal Infection 216 217 Society. Clin Orthop Relat Res. 2011; 469: 2992-2994.
- Workgroup Convened by the Musculoskeletal Infection Society. New definition for periprosthetic joint infection. J [18] 218 Arthroplasty. 2011; 26: 1136-1138. 219
  - Bedair H, Ting N, Jacovides C, Saxena A, Moric M, Parvizi J, Della Valle CJ. Diagnosis of early postoperative TKA infection using synovial fluid analysis. Clin Orthop Relat Res. 2011; 469: 34-40.
  - [20] Kim SG, Kim JG, Jang KM, Han SB, Lim HC, Bae JH. Diagnostic value of synovial white blood cell count and serum
  - C-reactive protein for acute periprosthetic joint infection after knee arthroplasty. J Arthroplasty. 2017; 32: 3724-3728. Yi PH, Cross MB, Moric M, Sporer SM, Berger RA, Della Valle CJ. Diagnosis of infection in the early postoperative period after total hip arthroplasty. Clin Orthop Relat Res. 2014; 472: 424-429.
  - Yu BZ, Fu J, Chei W, Hao LB, Chen JY. Neutrophil to lymphocyte ratio as a predictor for diagnosis of early periprosthetic joint infection. BMC Muscuoloskeletal Disorders. 2020; 21: 706.
  - Parvizi J, Gehrke T. International Consensus Group on Periprosthetic Joint Infection. Definition of periprosthetic joint infection. J Arthroplasty. 2014; 29: 1331.
- Sukhonthamarn K, Tan TI, Xu C, Ku FC, Lee MS, Citak M, Gehrke T, Goswani K, Parvizi J. Determining diagnostic [24] 230 thresholds for acute postoperative periprosthetic joint infection. J Bone Joint Surg Am. 2020; 102-A: 2043-2048. 231
- Fink B, Hoyka M, Weissbarth E, Schuster P, Berger I. The graphical representation of cell count representation: A new 232 procedure for the diagnosis of periprosthetic joint infections. Antibiotics (Basel). 2021; 10: 346. 233
  - Ghanem E, Parvizi J, Burnett RS, Sharkey PF, Keshavarzi N, Aggarwal A, Barrack RL. Cell count and differential of [26] aspirated fluid in the diagnosis of infection at the site of total knee arthroplasty. J Bone Joint Surg Am. 2008; 90-A:
- Spangehl MJ, Masri BA, O'Connell JX, Duncan CP. Prospective analysis of preoperative and intraoperative investigations 237 for the diagnosis of infection at the sites of two hundred and two total hip arthroplasties. J Bone Joint Surg Am. 1999; 238 239 81-A: 672-683.
  - McNally M, Sousa R, Wouthuyzen-Bakker M, Chen AF, Soriano A, Vogely HC, Lauss M, Higuera CA, Trebse R. The EBJIS definition of periprosthetic joint infection. A practical guide for clinicians. Bone Joint J. 2021; 103-B: 18-25.
- Choi H-R, Agrawal K, Bedair H. The diagnostic thresholds for synovial fluid analysis in late periprosthetic infection of 242 the hip depend on the duration of symptoms. Bone Joint J. 2016; 98-B: 1355-1359. 243
- Krenn V, Otto M, Morawietz L, Hopf T, Jakobs M, Klauser W, Schwantes B, Gehrke T. Histopathologic diagnostics in 244 endoprosthetics: periprosthetic neosynovialitis, hypersensitivity reaction, and arthrofibrosis. Orthopäde. 2009 Jun; 38(6): 245 520-30. 246
  - Krenn V, Morawietz L, Perino G, Kienapfel H, Ascherl R, Hassenpflug GJ, Thomsen M, Thomas P, Huber M, Kendoff D, Baumhoer D, Krukemeyer MG, Natu S, Boettner F, Zustin J, Kölbel B, Rüther W, Kretzer JP, Tiemann A, Trampuz A, Frommelt L, Tichilow R, Söder S, Müller S, Parvizi J, Illgner U, Gehrke T. Revised histopathological consensus classification of joint implant related pathology. Pathol Res Pract. 2014; 210: 779-786.
  - Müller M, Morawietz L, Hasart O, Strube P, Perka C, Tohtz S. Histopathological diagnosis of periprosthetic joint infection following total hip arthroplasty: Use of a standardized classification system of the periprosthetic interface membrane. Orthopade. 2009; 38: 1087-1096.
- Fink B, Hoyka M, Weissbarth E, Schuster P, Berger I. A new graphic type differentiation of cell account determination 254 255 for distinguishing acute periprosthetic joint infection from hemarthrosis. Antibiotics (Basel). 2022 Sep 21; 11(10): 1284.

12	B. Fink et al. / Graphic type differentiation of cell count data for diagnosis of early and late PJI
[34]	Tahta M, Simsek ME, Isik C, Akkaya M, Gursoy S, Bozkurt M. Does inflammatory joint diseases affect the accuracy of infection biomarkers in patients with periprosthetic joint infections? A prospective comparative reliability study. J Orthop Sci 2019; 24: 286-289.
[35]	Yi Ba PH, Corss MB, Moric M, Levine BR, Sporer SM, Paprosky WG, Jacobs JJ, Della Valle CJ. Do serologic and synovial tests help diagnose infection in revision hip arthroplasty with metal-on-metal bearings or corrosion? Clin Orthop Relat Res. 2015; 473: 498-505.
36]	Cooper HJ, Della Valle CJ, Berger RA, Tetreault M, Paprosky WG, Sporer S;, Jacobs JJ. Corrosion at the head-neck taper
37]	as a cause for adverse local tissue reactions after total hip arthroplasty. J Bone Joint Surg Am. 2012; 94-A: 1655-1661. Cooper HJ, Urban RM, Wixson RL, Meneghini RM, Jacobs JJ. Adverse local tissue reaction arising form corrosion at the femoral neck-body junction in a dual-taper stem with a cobalt-chromium modular neck. J Bone Joint Surg Am. 2013;
38]	95-A: 865-872. Earll MD, Earll PG, Rougeux RS. Wound drainage after metal-on-metal hip arthroplasty secondary to presumed delayed
39]	hypersensitivity reaction. J Arthroplasty. 2011; 26: 338.e5-7.  Mikhael MM, Hanssen AD, Sierra RJ. Failure of metal-on-metal total hip arthroplasty mimicking hip infection. A report
40]	Wyles CC, Larson DR, Houdek MT, Sierra RJ, Trsoudale RT. Utility of synovial fluid aspirations in failed metal-on-metal
	total hip arthroplasty. J Arthroplasty. 2013; 28: 818-823.
	40
	A Y
	XO .